


IDT Oligo ESI Analysis

Sales Order Number: 8281676

Reference ID: 96862228

Manufacturing ID: 178783626



Calculated Molecular Weight: 43,119.9

Measured Molecular Weight: 43,123.3

Sequence Name: Ultramer140-2

Oligo Sequence: 5'- CAC ATG ACC ATT TGG CGG TCT GTT CGC ACG ATA AAC
ACC GTG TAG AAT TAC GAT AAC TAT AAT ATC GCT CCC GGT
CTG TAA AAC TGC ATA TGG TGC TAC TCT AAG GGT ATG GGA
CAC TTC TTG AAT TCC TTG CAG TCT TG -3'